

Application recommendations:

Indication	Program
Degenerative and inflammatory diseases of the skeletal and motor apparatus (with the exception of acute episodes) Rheumatic diseases (arthrosis, arthritis, morbus Bechterew...), morbus Sudeck	P4
Inflammatory diseases of the skeletal and motor apparatus in acute episodes Rheumatic diseases	P1
Regulation of bone growth Osteoporosis, initial treatment of fractures subjected to conservative or surgical treatment, endoprosthesis loosening, pseudoarthroses, morbus Perthes	P2
Fibromyalgia	P4
Back pains cervical, thoracic and lumbar spine syndrome, ischias, lumbalgias,	P4
Sports injuries Acute joint and muscle traumas (muscle / ligament tears and strains, sprains, bruises, follow-up treatment of luxation, muscle affections, insertion tendinitis, epicondylitis, achillodynia,....	P4
Metabolic diseases, diseases of the digestive organs Ulcers of the stomach and duodenum, diabetes, gout	P4
Hypertension	P1
Hypotension	P3
Circulatory disturbances, lymph blockage Arterial occlusive disease, arteriosclerosis, venous stasis, haematomas, oedemas, angina pectoris, cardiac infarction,...	P3
Migraines Headaches, weather sensitivity	P1
Problems with falling asleep and sleeping through the night	P1
Wound healing	P4
Respiratory tract diseases Rhinitis, bronchial asthma, bronchitis, sinusitis, pulmonary emphysema, pneumonia,...	P1
Apoplexy	P3

Applications within the framework of general prevention and wellness:

Objective	Program
General old-age debility	P3
Calming, relaxation	P5
Exhaustion	P3
Enhancement of circulation	P3
Figure problems	P3
Climacteric complaints	P5
Muscle strains	P4
Menstrual pains	P5
Sleep disturbances	P1/P5
Disturbances of sexual function	P3
Stimulation, vigilance	P3

Contraindications

- ↳ Cardiac pacemakers, insulin pumps and other devices that are implanted or have to be worn next to the body constantly generally represent a contraindication for this magnetic field therapy. The therapeutic locus must be at a sufficient distance from the electrical or electronic implant.
Local treatments distal to the trunk are possible, but require the explicit permission of and stringent control by specialist medical staff.
- ↳ Pregnancy constitutes a contraindication. The reason for the contraindication is the insufficiency of the available scientific studies on treatment of pregnant women. The principle applied here is that of preventive risk avoidance.
- ↳ Small children should only be subjected to treatment with this device with the explicit permission of a physician.
- ↳ Treatment should be discontinued in acute conditions involving high fever.
- ↳ Degree III and IV cardiac dysrhythmias and convulsive disorders require monitoring by medical specialists.
- ↳ Patients who are allergic to nickel should not allow the "rod applicator" to contact their skin.

Important information

- ↳ The magnetic field has a range of about 1 metre. Therefore, persons in the immediate vicinity must also be taken into consideration.
- ↳ Metallic implants are not a problem with this therapy as they are in classic electrotherapy, since the pulsed magnetic field does not have a deep thermal effect.