

Program description

- 1 Calming program
- 2 Basic / Stimulation program
- 3 Vital and Fitness program
- 4 Regeneration/ Pain Relief program
- 5 Enzyme / Harmonization program

Program 1 - **8 Minutes** (recommended Therapy time) mostly around 4 Hz
Has a calming effect for Hypertension, Migraine, Weather sensitivity, Neurological problems, Nerve inflammation, Stomach / gastrointestinal problems, Ulcers, Heart problems and prophylaxis, Angina pectoris, Skin disease, faster Wound healing, burning sensation, Ulcus Cruris, Decubitus, pain of scars, Dermatitis, Neurodermatitis. All Acute inflammation and episodes in Rheumatic diseases. Insomnia. Respiratory tract diseases.

Program 2 - **16 Minutes** (recommended Therapy time) mostly around 20 Hz
Strengthening the all over improvement, and Immune system, Virus infection, strengthening the organism before and after chemotherapy.
Regulation of bone growth, Osteoporosis, bone fractures, endoprosthetic loosening, pseudoarthrosis.

Program 3 - **12 Minutes** (recommended Therapy time) mostly around 0,5-100 Hz
Vitalizing the organism. Chronique Fatigue (CFS-Syndrome), Exhausting sensations, Burn out Syndrome, a better well-being, Hypotension, Circulatory disturbances, Lymph-blockages, Arterial occlusive disease, Thrombosis and prevention, Arteriosclerosis, Haematoma, Oedemas, Lowering Follow up problems through immovability like Pneumonia und Thrombosis prophylaxis

Program 4 - **20 Minutes** (recommended Therapy time) oszcillating 2-22 Hz
Degenerative and inflammatory diseases of the Sceletal and motor apparatus(with exception of acute episodes), Arthrosis, Arthritis, Rheumatic diseases (M Bechterew, M Sudeck). Heart diseases, Heart insufficiency, Metabolic disturbances and diseases like Diabetes, Pancreas problems, Gout, Colitis etc. Fibromyalgia. Back pains.

Programm 5 - **20 Minutes** (recommended Therapytime) mostly between 0,5 - 1 Hz
Has been developed for stimulating the enzymes and is unique on the market. This Program should be used alone, or always in a combination with one of the above programs, where you can adjust the time of P5 on yourself.(4 min- 20 min)

Please be aware that you drink a large glass of water (0,5 liter) before and/or after the therapy.
This helps to regenerate faster and to remove the toxins out of the body.